

FM-02-0-014-1 Issue/Review Date: 1402.11.15	URS for Particle Counter	سینا دارو Sina Darou
Code No.: URS-QC-PC (1.0cfm)-001	Version No: Original	Approved Date: 1403.01.27
		Pages 1/7

User Requirement Specification (URS)






For
Particle Counter

Document Code: URS-QC-PC(1.0cfm)-001

☞ This document is under control and unalterable.

Revision History:

Change from previous version	Date
Original	1403.01.27

	Prepared By	Check By	Approved By	Approved By	Approved By
Organizational Title	Validation Supervisor	Laboratory Supervisor	QC Manager	QA Manager	Responsible Person
Name	E.Sedighi	A.Roostaie	Dr. Taghavi	Dr Keshavarz	Dr Torabi
Signature					
Date	1403.01.26	1403.01.26	1403.01.26	1403.01.27	1403.01.27




FM-02-0-014-1 Issue/Review Date: 1402.11.15	URS for Particle Counter	
Code No.: URS-QC-PC (1.0cfm)-001	Version No: Original	Approved Date: 1403.01.27
		Pages 2/7

Table of Content

1. Purpose and Scope	3
2. Responsibility	3
2.1 Preparation	3
2.2 Reviewing	3
2.3 Approval.....	3
2.4 Authorization	3
3. Overview	3
4. Operational Requirements	4
4.1 Function specification	4
4.2 Capacity	5
4.3 Process Requirements.....	5
4.4 Process Control	5
4.5 User Interface	5
4.6 Data and Security	5
5. Environment	5
6. Constraints	5
6.1 Milestones and Timelines	5
6.2 Equipment Constraints.....	5
6.3 Compatibility and Support.....	5
6.4 Availability/Maintenance	5
6.5 Access Level.....	5
7. Lifecycle	6
7.1 Development.....	6
7.2 Testing	6
7.3 Delivery/ Documentation.....	6
7.4 Support and training	7
8. Support and training	7
9. References	7
10. Definitions, Acronyms and Abbreviations.....	7
11. Document Distribution	7



FM-02-0-014-1 Issue/Review Date: 1402.11.15	URS for Particle Counter	
Code No.: URS-QC-PC (1.0cfm)-001	Version No: Original	Approved Date: 1403.01.27
		Pages 3/7

1. Purpose and Scope

This standard establishes a test procedure for evaluating the performance of air cleaning devices as a function of particle size. This document is applicable on testing for spot-checking real-time cleanliness during operation, root out contamination and creating standardized testing routine for long-term monitoring and data logging.

2. Responsibility

1. The expert of IPQC unit is responsible for preparing URS.
2. The manager and head of the QC are responsible for the approval of this document.
3. The quality assurance manager and technical officer are responsible for approving this document.

3. Overview

Optical particle counters (OPC) or spectrometers employ a small sensing volume, either by a focused incandescent lamp or by a laser source. In such instruments it is important to avoid coincidence errors resulting from more than one particle in the sensing volume. The instrument manufacturer specifies the maximum number concentration which can be handled. Generally, commercial instruments handle concentrations of up to $\sim 10^6$ particles/litre. Beyond this concentration limit, sample dilution is usually used, which decreases the accuracy of the determination of the concentration. The range of particle diameter that single-particle instruments are capable of handling is $\sim 0.1 \mu\text{m}$. Optical particle counters have found wide use, first in cleanroom monitoring and more recently in community air pollution and industrial hygiene studies. A number of laboratory instruments employing single-particle scattering have been constructed. A critical review of such instruments is given by Chigier and Stewart.

Based on previous light scattering experiments, the development of airborne particle counters and sizers started in the middle of the twentieth century and accelerated in the 1960s after the invention of lasers.

The particle size distribution in single-particle counters is determined by comparing the detected pulse heights of the optical signals that correspond to the single particle flow through a small illuminated zone with a standard calibration curve, obtained from a set of uniform particles of known diameter. For determination of the concentration, the flow rate (the volume to be tested) is also measured simultaneously. The different light detection geometries applied in different particle counters yield a wide range of instrument designs and constructions. As the dependence of the scattered intensity on the refractive index of the particles is less pronounced around 90° of scattering angle, in most airborne particle counters perpendicular scattering geometry is implemented.

FM-02-0-014-1

Issue/Review Date:
1402.11.15

URS for

Particle Counterسینا دارو
Sino Darou

Code No.:

URS-QC-PC (1.0cfm)-001

Version No:

Original

Approved Date:

1403.01.27


Pages

4/7

4. Operational Requirements**4.1 Function specification**

Size Channels	0.3 μm , 0.5 μm , 1.0 μm , 3.0 μm , 5.0 μm , 10.0 μm
Zero Count	1 count or less in 5 minutes
Coincidence Loss	Less than 5% at 400000 particles/ ft^3
Flow Rate	1.0 CFM
Counting Efficiency	50 % for 0.3 μm
Light Source	Long Life Laser™ diode with 10-year Mean Time to Failure (MTTF)
Count display	LCD
Maximum Count	9999999 display
Pump Type	Air vacuum, rated for continuous use
Sample /Hold Times	1 second to 24 hours
Count Alarm	1 to 9999999 counts
Count Cycles	Up to 100 while in automatic mode
Location label	0 to 999, appears on printout
Data Storage	5,000 samples, scrollable on Historical Data review screen
Battery	Rechargeable Li
Operating Time (Battery)	up to 6 hours
Battery Recharge Time	7 hours minimum
Communication	RS-232, RS-485, 802.11 b/g Wireless, 802.3 Ethernet, USB Client, USB Host
Weight	Max 7.0 Kg




FM-02-0-014-1 Issue/Review Date: 1402.11.15	URS for <h2 style="text-align: center;">Particle Counter</h2>	
--	--	---

Code No.: URS-QC-PC (1.0cfm)-001	Version No: Original	Approved Date: 1403.01.27	Pages 5/7
--	--------------------------------	-------------------------------------	---------------------

4.2 Capacity N.A
4.3 Process Requirements N.A
4.4 Process Control N.A
4.5 User Interface N.A
4.6 Data and Security N.A
5. Environment This device is placed in cleanroom where you'll find particle counters.
6. Constraints
6.1 Milestones and Timelines N.A
6.2 Equipment Constraints N.A
6.3 Compatibility and Support N.A
6.4 Availability/Maintenance N.A
6.5 Access Level N.A



FM-02-0-014-1 Issue/Review Date: 1402.11.15	URS for Particle Counter	
Code No.: URS-QC-PC (1.0cfm)-001	Version No: Original	Approved Date: 1403.01.27
		Pages 6/7

6.6 Procedural Constraints

- Since the handheld particle counter is small, it does not force airflow. Thus, due to the lower airflow, it can take longer than other particle counter options to produce reading.
- A handheld particle counter is designed for spot checks, so it does not integrate well into a real time system.

7. Lifecycle

7.1 Development


.N.A

7.2 Testing

- Verify that all aspects of the cleanroom system which contribute to its operational integrity (air handling, filtration systems, walls, ceilings, floors, etc.) are complete and functioning. The following primary tests shall be completed prior to performing the cleanliness classification tests:
 - Airflow test
 - Filter leak test
 - Room pressurization test
- The Cleanroom has been purged for at least 12 hours before testing commences.
- The rooms are divided into test grids as per minimum sampling locations requirement by ISO 14644-1:2015.
- The counter's sample rate is set to 1 m³/min and the sampling time set to 1 minute or as required.
- A minimum of 1 sample is taken at each grid, at a height of 1 meter above the floor.
- When an obstruction is encountered, the sample is to be taken at 300 mm above the obstruction.

7.3 Delivery/ Documentation

The manufacturer must comply with good documentation practice and deliver the necessary documentation to prove that the facility, system and equipment is constructed according to GMP

FM-02-0-014-1	URS for	Particle Counter		
Issue/Review Date: 1402.11.15				
Code No.: URS-QC-PC (1.0cfm)-001	Version No: Original	Approved Date: 1403.01.27	Pages 7/7	

guidelines and in line with our URS. Thus, it must be documented that the materials used are compliant with GMP guidelines and the appropriate validation tests have been performed. Checklists should be used to validate that all parts of construction, equipment and systems are in concordance with GMP guidelines. These should include, but not be limited to: the material used, quality specification, certificates, test methods/protocols, responsible test person, number of tests performed, test results, transfer plans and reports, IQ, OQ, PQ, Certificate of Analyses (CoA)/Compliance (CoC), validation and documentation

7.4 Support and training

Training course should be performed by manufacturer on the site.

8. References

- ISO 14644-1:2015(E) Cleanrooms and associated controlled environments-Part 1: Classification of Air Cleanliness by Particle Concentration.
- ISO 21501-4:2018 Determination of particle size distribution — Single particle light interaction methods — Part 4: Light scattering airborne particle counter for clean spaces

9. Attachments

N.A

10. Definitions, Acronyms and Abbreviations

N.A

11. Document Distribution

Department	Version
Original version	QA
Copy No. 1	QC

